

B-Life

Silicone face mask Kit

6 sizes



Autoclavable face masks suitable for reanimation and assisted ventilation.

Specific features

Connection standard ISO 15/22 mm
 Transparent shell improve patient monitoring
 Soft and hypoallergenic silicon pillow
 Latex Free

Technical data

Materials	Silicone, PC
Max autoclave temperature	134°C
Weight	300 g

Masks included

RM20400B	B-Life silicone mask sz.0
RM20401B	B-Life silicone mask sz.1
RM20402B	B-Life silicone mask sz.2
RM20403B	B-Life silicone mask sz.3
RM20404B	B-Life silicone mask sz.4
RM20405B	B-Life silicone mask sz.5

Class IIa Med.Device according to Dir. 93/42/CEE

B-LIFE KIT SIX SIZES –
 SILICONE FACE MASKS



RM20410B

CND Classification R03010199

Registration number 145641

Optionals

BG02001A	Laxan transparent bag for manual resuscitators
----------	--

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it