

Eco B-Life

Face mask



Disposable face mask suitable for reanimation and assisted ventilation.

Specific features

Transparent shell improve patient monitoring
Soft and hypoallergenic silicon pillow

Technical data

Materials	Rubber, PC
Dimensions	110x90x70 mm
Weight	61 g

Optionals

BG02001A	Laxan transparent bag for manual resuscitators
----------	--

Class IIa Med.Device according to Dir. 93/42/CEE

ECO B-LIFE SIZE 3 - RUBBER
FACE MASK



RM20510B

CND Classification R03010199

Registration number 104840

Rev.0 (18/04/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
spencer.it