

D-COLLIN

Tongue forceps



Traditional stainless steel tongue forceps.

Specific features

The gear rack allows multiple selection of closure pressure
Satin finish

Technical data

Length	160 ± 5 mm
Weight	48 ± 5 g
Gear rack positions	6
Material	Stainless steel

Class I Med.Device according to Regulation 2017/745/UE

D-COLLIN - STAINLESS STEEL TONGUE
FORCEPS

RM30301

CND Classification



V9099

Distributed by Spencer Italia s.r.l.

Rev.0 (25/08/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it