

RVX 02

Suction canister for Res-Q-Vac



Autoclavable canister to be used as spare part for manual suction device Res-Q-vac.

Specific features

Integrated overflow valve
High transparency
Sterilizable at 121°C/1 bar

Technical data

Dimension	ø92x115 mm
Capacity	380 ml
Material	PC
Weight	160 g

Class I med.dev. according to Dir. 93/42/CEE

SUCTION VESSEL AUTOCLAVABLE
LT.0.380



RV71102A

CND classification	V9099
--------------------	-------

Registration identifier	818169
-------------------------	--------

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it