

VP FILTER



Disposable filter designed to protect patient and suction device from cross infections. Prevents passage of fluids that come in contact with it.

Specific Characteristics

The filter is hydrophobic and provides both viral and bacterial filtration

Easy to install

The device is unidirectional; input and output ports are clearly marked

Technical data

Dimensions	ø65x 50 mm
Weight	20 g
Input and output Connection	Hose ø 7 mm <small>(internal 3,7 mm)</small>

VP FILTER - FILTER FOR ASPIRATORS
RV71106C

Distributed by Spencer Italia S.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it