## **B-VAC**

Manual aspirator for secretions with adult catheter





Manual aspirator for secretions with reduced obstruction capacity, its function is the removal of fluids from respiratory tracts or from the respiration support system and of infective material from wounds. The aspirator creates a depression which pulls fluids out through the disposable catheter tied to the disposable vessel. It does not need power source

## Specific features

- Lightweight and compact
- Usable with only one hand
- Disposable vessel with integrated filter Selector for adult or paediatric use

## Accessories/Consumed replacement parts

- RV71151E KIT disposable vessel c/adult tube
- RV71156E paediatric suction tube
- RV71157E adult suction tube

## **Compliant with**

EN ISO 10079-2

we should also a	
Technical data	
Dimensions with vessel	210x170x75 mm
Adult catheter dimensions	ø16x22 cm
Max Vacuum	550 mmHg
Flow	>20 l/min
Vessel's capacity	350 ml
Viral and bacterial filtration efficiency	99,99%
Materials	PE, PVC
Aspirator's weight w/vessel	270 g
Weight adult catheter	30 g
Vacuum selector	100% o 50%
Class I Med. Device according to Dir. 93/42/CEE	
B-VAC - MANUAL ASPIRATOR FOR	
SECRETIONS	
RV71150B	
CND Classification	V9099

**Registration number** 

277342

Rev.2 (15/01/2018) UNCHECKED COPY – further revisions will be available on <u>http://support.spencer.it</u>

Spencer Italia S.r.I. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it