

## D-B-VAC

### Manual aspirator for secretions with adult catheter



Manual aspirator for secretions with reduced obstruction capacity, its function is the removal of fluids from respiratory tracts or from the respiration support system and of infective material from wounds. The aspirator creates a depression which pulls fluids out through the disposable catheter tied to the disposable vessel. It does not need power source

#### Specific features

- Lightweight and compact
- Usable with only one hand
- Disposable vessel with integrated filter
- Selector for adult or paediatric use

#### Accessories/Consumed replacement parts

RIRV005	Spare disposable canister
RIRV006	Paediatric suction tube
RIRV007	Adult suction tube

#### Compliant with

EN ISO 10079-2

#### Technical data

Dimensions with vessel	210x170x75 mm
Adult catheter dimensions	ø16x22 cm
Max Vacuum	550 mmHg
Flow	>20 l/min
Vessel's capacity	350 ml
Viral and bacterial filtration efficiency	99,99%
Materials	PE, PVC
Aspirator's weight w/vessel	270 g
Weight adult catheter	30 g
Vacuum selector	100% o 50%

Class I Med.Device according to Regulation 2017/745/UE

B-VAC - MANUAL SUCTION DEVICE

RV71170B



CND Classification

V9099

Distributed by Spencer Italia s.r.l.

Rev.0 (24/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)