

BLANCO

Portable suction device with battery



Portable surgical Aspirator designed for the aspiration of bodily fluids in adults and children. It is particularly suitable for nasal, oral and tracheal aspiration of mucus, catarrh or blood after minor surgical procedures and can be used in post-operative therapy at home or conveniently transported from one hospital ward to another.

Specific features

- Compact and lightweight
- Fast reading of the vacuum level thanks to the dimensions of the vacuum gauge
- Classification ISO High vacuum/Low flow
- Comfortable transport handle
- ABS body easy to clean
- Led indicator for battery status
- Multiple power supplies make it extremely versatile

Standard equipment

- 1000 ml autoclavable canister
- Silicon suction tubes
- Antibacterial filter
- Sterile suction probe
- Power adapter
- 12V car adapter

Accessories/Spare parts

- | | |
|----------|-----------------------------|
| SC70085A | Autoclavable canister 1l |
| SC70088A | Disposable conical junction |
| SC10037C | Antibacterial filter |
| SC10014C | Carry bag for Blanco |
| SC10012B | 10G Fastener for Blanco |

Compliant with

- ISO 10079-1
- EN 60601-1
- EN 60601-1-2
- EN 60601-1-11

Technical data

Dimensions	350x210x180 mm
Weight	3,5 kg
Max Vacuum ^(without canister)	0,75 bar
Max Flow ^(with canister)	16 l/min
Operating mode	Continuous
Power supply	230V/50hz
	P.Adapter 100-240V 50-60hz
	Internal battery (12V 4A) Car adapter (12V 4 A)
IP Protection	IP 21
Insulation class (When operated through adapter)	Class II
Battery operating time	80 min
Battery charging time	240 min
Noise level	63 db
Operating temperature	5°C / 35°C
Storage temperature	-25°C / 70°C

Class IIa Med.Device according to Dir. 93/42/CEE

BLANCO TRANSP. 12V 230 16 L/MIN
W/BATTERY
SC10002B



CND Classification	Z120105
Registration number	307630

Rev.3 (27/06/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it