

# BLANCO

## Portable suction device



Suction device operating at 12V, designed for the aspiration of bodily fluids in adult and children. It's particularly suitable for nasal, oral or tracheal aspiration of mucus, catarrh or blood after minor surgical procedures and can be used in post-operative therapy at home or conveniently transported from one hospital ward to another.

The device is intended to be powered by a 12V DC car battery with non continuous operation (Ton: 20 min / Toff: 40 min).

### Specific features

- Compact and with low weight
- Fast reading of the vacuum level thanks to the dimensions of the vacuum gauge
- Classification ISO High vacuum/High flow
- Comfortable transport handle
- ABS body easy to clean

### Standard Equipment

- 1000 ml autoclavable canister
- Silicon suction tubes
- Antibacterial filter
- Sterile suction catheter
- Cigarette lighter adapter

### Accessories/Spare parts

- |          |                                      |
|----------|--------------------------------------|
| SC70085A | Autoclavable canister 1l             |
| SC70088A | SCX 88 – Disposable conical junction |
| SC10037C | Antibacterial filter                 |
| SC10014C | Carry bag for Blanco                 |
| SC10012B | 10G Fixing for Blanco                |

### Compliant with

- ISO 10079-1
- EN 60601-1
- EN 60601-1-2

### Technical data

- |                                  |                       |
|----------------------------------|-----------------------|
| Dimensions                       | 350x210x180 mm        |
| Weight                           | 2,5 kg                |
| Max Vacuum (without canister)    | 0,75 bar              |
| Max Flow (without canister)      | 25 l/min              |
| Operating mode                   | 20 min ON/ 40 min OFF |
| Power supply                     | 12V cc                |
| Power consumption                | 40W                   |
| Electrical safety classification | Class II type B       |
| Operating temperature            | 5°C / 35°C            |
| Storage temperature              | -25°C / 70°C          |

Class IIa Med.Device according to Dir. 93/42/CEE

BLANCO 12 V 25 L/MIN -  
PORTABLE SUCTION DEVICE

SC10003B



- |                     |         |
|---------------------|---------|
| Classification CND  | Z120105 |
| Registration number | 307643  |

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)