

BLANCO

Suction device



Blanco is a desk-type electric suction unit for aspiration of body liquids, nasal, oral and tracheal aspiration in adults or children. Particularly suitable for its particularly suitable for nasal, oral or tracheal aspiration of mucus, catarrh or blood after minor surgical procedures and can be used in post-operative therapy at home or conveniently transported from one hospital ward to another. Easily portable equipment designed for continuous use.

Specific features

- Compact and lightweight
- Fast reading of the vacuum level thanks to the dimensions of the vacuum gauge
- Classification ISO High vacuum/low flow
- Comfortable transport handle
- ABS body easy to clean

Standard Equipment

- Vaso autoclavabile da 1000 ml
- Silicon suction tubes
- Antibacterial filter
- Sterile suction catheter
- Shuko power cable

Compliant with

- ISO 10079-1
- EN 60601-1
- EN 60601-1-2

Dati tecnici

Dimensions	350x210x180 mm
Weight	3,6 kg
Max Vacuum <small>(without canister)</small>	0,80 bar
Max Flow <small>(without canister)</small>	40 l/min
Operating mode	Continuous
Power supply	230V/50hz
Power consumption	110 VA
IP protection	IP 21
Insulation class	Classe II tipo BF
Operating temperature	5°C / 35°C
Storage temperature	-25°C / 70°C

Class IIa Med.Device according to Dir. 93/42/CEE

BLANCO – PORTABLE SUCTION
DEVICE 220V 40L/min
SC10004B



CND Classification	Z120105
Registration number	307650

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it