

BLANCO FASTENER

10 G Fastener for Blanco suction units



Blanco fastener fixing systems are intended for installation in ambulances to anchor Spencer suction devices.

Specific features

- Stainless steel structure confers high resistance
- The fixing is done simply by turning the knob
- The smooth surfaces are easy to clean

Technical data

Dimensions (DxLxH)	258 x 250 x 232 ± 5 mm
Materials	Stainless steel, Nylon
Weight	1.10 ± 0.03 kg

Class I MD compliant with UE Reg. 2017/745

FIXATION FOR BLANCO SUCTION
DEVICE 10G
SC10012B



Rev.0 (30/11/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it