

## Blanco fastener



Fastener to allow safe wall installation of Blanco suction devices.

### Specific Characteristics

Stainless steel structure confer high resistance  
 The fixing is done simply by turning the knob  
 The smooth surfaces are easy to clean

### Technical data

Dimensions	258x250x232 ± 5mm
Materials	Steel, Nylon
Weight	1,10 ± 0,03 kg

Class I Med.Device according to Dir. 93/42/CEE

FIXATION FOR BLANCO SUCTION  
 DEVICE  
**SC10013B**



Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)