

# BLANCO

## Portable suction device



Blanco is a desk-type electric suction unit for aspiration of body liquids, nasal, oral and tracheal aspiration in adults or children. The appliance is designed for easy transport and almost continuous use. The unit is operated by means of an electronic power management system, also monitoring the internal battery operation and status. A LED system placed on the front panel indicates appliance start-up and battery charge status. The unit is equipped with suction regulator on the front panel and polycarbonate autoclavable jar complete with overflow valve. Thanks to these characteristics and to its functions, this device is particularly suitable for different applications: utilization in hospital wards, for tracheotomy, for suction of body liquids and for minor surgery.

### Specific Characteristics

Compact and lightweight
Fast reading of the vacuum level thanks to the dimensions of the vacuum gauge
Classification ISO High vacuum/low flow
Comfortable transport handle
ABS body easy to clean
Not replaceable internal rechargeable battery
Led indicator for battery status

### Standard Equipment

1000 ml autoclavable canister
Silicon suction tubes
Antibacterial filter
AC/DC universal adapter with detachable power cord
12V car adapter

### Accessories/Spare parts

SC70085A	Autoclavable vase 1l
SC70088A	SCX 88 – Disposable conical junction
SC10037C	Antibacterial filter
SC10014C	Carry bag for Blanco
SC10012B	10G Fixing for Blanco

### Compliant with

ISO 10079-1
EN 60601-1
EN 60601-1-2

### Technical data

Dimensions	350x210x180 mm
Weight	4,39 kg
Max Vacuum <sup>(without canister)</sup>	0,80 bar
Max Flow <sup>(without canister)</sup>	36 l/min
Max Flow <sup>(with canister)</sup>	23 ± 2 l/min
Operating mode	Continuous
Power supply	14Vcc 4A with AC/DC adapter (input: 100-240 V~ - 50/60 Hz - 100 VA) – Internal rechargeable Pb battery (12Vcc 4A) – Car adapter (12V – 4A)
Power consumption	4 A
IP protection	IP 21
Insulation class	Class I/ Class II / or internally powered
Applied part	Type BF
Noise level	Max 68,5 db
Battery operating time	60 min
Battery charging time	120/150 min
Operating temperature	5°C / 35°C
Storage temperature	-40°C / 70°C

### Class IIa Med.Device according to Dir. 93/42/CEE

BLANCO PORTABLE ASPIRATOR 12V 36 L/min W/BATTERY  
SC10036B



Classification CND	Z120105
Registration number	841047

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)