

BLANCO

Portable suction device



Blanco is a desk-type electric suction unit for aspiration of body liquids, nasal, oral and tracheal aspiration in adults or children. The appliance is designed for easy transport and almost continuous use. The unit is operated by means of an electronic power management system, also monitoring the internal battery operation and status. A LED system placed on the front panel indicates appliance start-up and battery charge status. The unit is equipped with suction regulator on the front panel and autoclavable jar complete with disposable bag. Thanks to these characteristics and to its functions, this device is particularly suitable for different applications: utilization in hospital wards, for tracheotomy, for suction of body liquids and for minor surgery.

Specific Characteristics

Compact and lightweight
Fast reading of the vacuum level thanks to the dimensions of the vacuum gauge
Classification ISO High vacuum/low flow
Comfortable transport handle
ABS body easy to clean
Not replaceable internal rechargeable battery
Led indicator for battery status

Standard Equipment

1000 ml autoclavable canister and disposable bag
Silicon suction tubes
AC/DC universal adapter with detachable power cord
12V car adapter

Accessories/Spare parts

SC75017D	Oval canister for Serres Bag
SC75018E	36 Serres bags
SC10014C	Carry bag for Blanco
SC10012B	10G Fixing for Blanco

Compliant with

ISO 10079-1
EN 60601-1
EN 60601-1-2

Technical data

Dimensions	350x245x180 mm
Weight	4,45 kg
Max Vacuum <small>(without canister)</small>	0,80 bar
Max Flow <small>(without canister)</small>	36 l/min
Operating mode	Continuous
Power supply	14Vcc 4A with AC/DC adapter (input: 100-240 V~ - 50/60 Hz - 100 VA) – Internal rechargeable Pb battery (12Vcc 4A) – Car adapter (12V – 4A)
Power consumption	4 A
IP protection	IP 21
Insulation class	Class I/ Class II / or internally powered
Applied part	Type BF
Noise level	Max 68,5 db
Battery operating time	60 min
Battery charging time	120/150 min
Operating temperature	5°C / 35°C
Storage temperature	-40°C / 70°C

Med. Dev. CE marked according to Dir. 93/42/CEE

MOD.D018 BLANCO – PORTABLE
ASPIRATOR 36L W/DISPOSABLE BAG
SC10039B



Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it