

## SCX 85

Autoclavable canister for suction devices



Autoclavable canister with overflow valve to be used with suction devices.

### Specific features

- Autoclavable up to 30 cycles
- The overflow valve prevents spillage of fluids if the maximum capacity is reached
- High transparency to allow proper view of the content

### Technical data

Canister capacity	1000 ml
Weight	330 g
Material	PC
Autoclave settings	15min at T 121°/ P.rel 1 bar

### Related products

SC75016E	Sterile patient suction tube 1,8 m
OX10105A	S-TUBE – Silicone tube ø8x14mm
SC10037C	Spare filter for Blanco
RISC005E	Disposable filter for Ambujet - Homejet

Class I med. Dev. according to dir.93/42/CE  
 SUCTION VESSEL AUTOCLAVABLE LT.1  
 W/ SPILL OVER PRO  
**SC70085A**



Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)