

SCX 88

Conical fitting



Disposable conical fitting to be used with Spencer Blanco suction devices.

Specific Characteristics

The shape allow the use with suction tubes of different diameters

Easy to apply

Suitable for silicon tubes with internal diameter 8,9 and 10 mm

Technical data

Dimensions	Ø15,3x58,5 mm
Internal diameter	4,9 mm
Weight	4 g

SCX 88 - CONICAL JUNCTION

SC70088A

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it