

ANTIBACTERIAL FILTER



Disposable antibacterial filter for Spencer suction devices. Its application is fundamental to prevent the passage of contaminants between the collection canister and the suction device.

Specific Characteristics

The transparent body allows to see if any fluids are inside the filter

The flow direction is clearly marked on the filter, allowing a proper application

Technical data

Dimensions	68x41x71 mm
Weight	13 g
Connection	ø 10,5/13,8 mm

DISPOSABLE FILTER FOR VACU AIDE 2000 AND JET
SC73013E

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it