

## Sterile connecting tube



Disposable connecting tube with female conical end parts, intended for use with suction devices.

### Specific Characteristics

Transparent to allow to see the fluid passage  
Flexible conical female connectors  
Provided in sterile package – Sterilized by EO  
Contains DEHP Phtalates

### Technical data

Internal diameter	7,2 mm
External diameter	9,8 mm
Length	1,8 m
Material	PVC
Weight	120 g

Class IIa Med.Device according to Dir. 93/42/CEE

PATIENT SUCTION TUBE 1,8 MT FOR  
VACU AIDE/JET  
SC73016E



Distributed by Spencer Italia S.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)