

## Reusable collection canister



Reusable collection canister with overflow protection.

### Specific features

- Sterilizable up to 30 times
- The overflow protection prevents suction when the maximum level is reached
- High transparency to allow view of the content

### Technical data

Canister capacity	1000 ml
Weight	330 g
Material	PC
Sterilization parameters	15min at T 121° / Press.1 bar

### Related products

SC75016E	Sterile suction tube 1,8 m
OX10105A	S-TUBE – SILICONE SUCTION TUBE. ø8x14mm
SC10037C	Disposable filter for Blanco
RISC005E	Disposable filter for Ambujet - Homejet

Class I Med.Dev. according to dir.93/42/CE

1LT COMPLETE REUSABLE CANISTER T

SC73067A



Rev.0 (14/05/2019)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)