

Reusable canister



Reusable collection canister with overflow protection.

Specific features

Sterilizable up to 30 times at 121°C
 The overflow protection prevents suction when the maximum level is reached
 Silicon seal ring

Technical data

Canister capacity	1000 ml
Material	PC/ABS

Class I Med.Dev. according to dir.93/42/CE
 1000 ML REUSABLE CANISTER W/LID
 SC73069A



Distributed by Spencer Italia s.r.l.

Rev.0 (23/07/2021)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it