

300 ml Spare canister

For Jet Compact



Disposable canister with suction tube for use with suction device Jet Compact..

Specific features

- Integrated hydrophobic filter
- Supplied with disposable suction tube
- High transparency to allow proper view of the content

Technical data

Canister capacity	300 ml
-------------------	--------

Class I med. Dev. according to dir.93/42/CE
 DISPOSABLE JAR 300 ml FOR JET WITH
 PATIENT TUBE
 SC75001E



Distributed by Spencer Italia s.r.l.

Rev.0 (26/07/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it