

BATTERY FOR JET



Rechargeable battery intended for use with Spencer Jet series suction devices.

Specific Characteristics

Lightweight and compact

The shape of the plug ensures a correct and safe connection

Compatible with the suction devices HomeJet, AmbuJet, Jet Wire, Jet Compact

Technical data

Dimensions	^{Padding excluded} 72x30x54 mm
Weight	280 g
Type	NiMH
Voltage	12 V
Capacity	1,6 Ah
Protection systems	Fuse and thermistor

12V Ni-MH BATTERY FOR JET
SC75007K

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it