

Carry Bag for Jet Compact



Carry bag to be used with the suction device Jet Compact.

Specific Characteristics

- PVC Fabric easy to clean
- The top part transparency, allow an easy visualization of the setted vacuum
- The replacement of the canister and battery charging operations, can be made without removing the device from the bag
- Elastic to hold the canister safely
- Anterior pocket for housing of suction tubes
- Shoulder belt with aluminum carabiners to make the transport easier

Technical data

Dimensions with Jet inside	80x190x190 ± 5mm
Materials	PVC, PP
Shoulder belt length	50 cm
Shoulder belt carabiners	Anodized Al
Overall weight	230 g

TRANSPORT BAG FOR JET COMPACT
SC75011E

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it