

## SERRES CANISTER

With disposable bag



Collection jar with disposable bag for collection of secretions. The compatibility with some aspirators can be obtained with the use of specific conversion kits.

### Specific Features

- Steam sterilizable canister
- Disposable collection bags with integrated antibacterial filter
- The antibacterial filter acts also as an overflow valve, preventing spillage when the maximum capacity is reached
- High transparency to allow easy view of the content

### Technical data

Canister capacity	1000 ± 100 ml
Weight w/out bag	270 g
Weight of the bag	50 g

### Spare parts

SC75017D	Serres canister
SC75018E	36 disposable bags

Canister and bag CE marked according to dir.93/42/CE  
 SERRES CANISTER 1000 ml W/BAG  
**SC75015D**

Distributed by Spencer Italia S.r.l.

Rev.1 (15/01/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)