

## SERRES CANISTER



Collection canister for suction devices to be used with dedicated disposable bags. Compatibility with some aspirators can be achieved with the use of dedicated conversion kits.

### Specific features

Steam sterilisable  
High transparency to allow easy view of the content

### Technical data

Canister capacity	1000 ± 100 ml
Weight	270 g

### Accessories

SC75018E	BOX 36 PCS SERRES BAGS 1000 ml
----------	--------------------------------

Canister CE marked according to dir.93/42/CE

1000ml OVAL CONTAINER FOR SERRES  
BAGS

**SC75017D**

Distributed by Spencer Italia S.r.l.

Rev.0 (09/10/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)