

Disposable bags



Package with 36 disposable bags intended to be used with Serres collection canisters.

Specific Characteristics

Transparent to allow proper view of the content
Easy to install

Technical data

Capacity	1000 ± 100 ml
Weight of the bag	50 g

Related devices

SC75017D	Collection canister Serres 1000 ml
----------	------------------------------------

Class I Med.Device according to Dir. 93/42/CEE

PACKAGE 36pcs 1000ml SERRES
DISPOSABLE BAGS
SC75018E



Classificazione CND

A060399

Distributed by Spencer Italia S.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it