

Jet conversion Kit



Conversion Kit for Home Jet and Ambu Jet suction devices intended to allow the housing of Serres type canisters.

Specific Characteristics

Autoclavable canister
Disposable suction bags
High transparency to allow adequate visualization of the content
Easy to install

Technical data

Canister capacity	1000 ± 100 ml
Canister weight	270 g
Single disposable bag weight	50 g
Weight of the kit	2,2 kg

Spare parts

SC75017D	Serres canister with connector
SC75018E	Pack of 36 disposable bags

Canister and bags are CE marked according to dir.93/42/CE
 CONVERSION KIT AMBUJET/HOMEJET
 W/SERRES
SC75021D

Canister and bags are distributed by Spencer Italia S.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it