AMBUJET P 800 D 10G

Portable suction device with 10G wall fastener with auto recharge





Space-saving suction device to be used to remove fluids from the airway or respiratory support system. The device creates a negative pressure that draws fluids through disposable tubing that is connected to the collection container. The base performs as a fastener in accordance to EN 1789 standard.

Specific Characteristics

- Low weight and size
- Easy view of the vacuum set thanks to the led indicators Classification ISO High Vacuum/ High Flow
- Suction device protected from bumps and falls thanks to the protective shell
- Mechanism of engagement and release from the base extremely rapid
- Direct connection between canister and vacuum pump Recharging of the battery is done simply by placing the device on the mounting system

Standard Equipment

- 10g fastener/charger with quick unlocking system Battery
- Disposable canister
- Disposable patient tube 1,8 m
- Shoulder strap for easy transport

Optional accessories/Consumables

SC75012E	Lateral bags for Jet
SC75014D	Supply cable for cig.lighter 12vcc
SC73017E	220V Power supply
RISC005E	Antibacterial filter
SC73016E	Sterile suction catheter 1,8m
SC73039E	800 ml disposable canisters (10pcs set)

Compliant with

ISO 10079-1 EN 1789

Technical data			
Jet Dimensions	200x225x285 mm		
Jet Dimensions w/fastener	200x250x320 mm		
Suction device weight	2,4 ± 0,1 kg		
Fastener weight	0,8 ± 0,1 kg		
Vacuum Range	From 50 to 550 mmHg		
Flow (free flow)	32 l/min		
Noise level	70 dB		
Electrical connection requirements	100-240 VAC 47-63 Hz 0,75A max; 12VDC 33W max		
IP Grade	IP12		
Battery	12V DC Ni-MH		
Operating time	30 min. ON, 30 min. OFF		
Battery recharging time	5 h		
Container capacity	800 ml		
Patient tube lenght	1,8 m		
Operating temp.	From 0°C to +45°C		
Storage temp.	From – 40°C to +70°C		
Umidity	From 0 to 95%		

Med.Device CE marked according to Dir. 93/42/CEE			
AMBUJET P 800 D 10G			
SC75550C			
Registration number	81583		

Rev.0 (31/07/2018) UNCHECKED COPY – further revisions will be available on http://support.spencer.it Spencer Italia S.r.I. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it