


PILLOW FOR SUPER BLUE



Spare pillow for head immobilizer Super Blue.

Specific Characteristics

Wide strap surfaces make the headblock positioning easier

Pillow made of non-absorbent material for an easy sanitation

Large holes for ear inspection

Technical data

Pillow dimensions	255x85x160 mm
Hole for ear inspection	ø70 mm
Pillow weight	500 ± 10g
Material	PE

CUSHION FOR SUPER BLUE WITH
SPENCER LOGO
SH00108C

Rev.1 (15/03/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it