

SUPER

Compact universal head immobilizer



Super head restraints are accessories for spine boards or pick-up stretchers, to be used to increase the degree of patient immobilisation.

Specific features

- The padding of the base increase patient comfort
- Wide strap surfaces make the headblock positioning easier
- Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
- Pillows made of non-absorbent material for an easy sanitation
- Large holes for ear inspection
- Allows for preliminary diagnostic examinations with X-rays

Technical data

| | |
|--|-------------------|
| Base dimensions ^(buttonhole excluded) | 365 x 260 x 20 mm |
| Pillow dimensions ^(1 piece) | 255 x 85 x 160 mm |
| Hole for ear inspection | Ø 70 mm |
| Base weight | 360 ± 10 g |
| Pillow weight | 500 ± 10 g |
| Overall weight | 1,46 ± 0,07 kg |
| Materials | PE, PVC, Nylon |

Class I MD compliant with UE Reg. 2017/745

SUPER BLUE BLUE - COMPACT
UNIV. HEAD IMMOB.
SH00110C



| | |
|---------------------|---------|
| CND Classification | Y060312 |
| Registration number | 102238 |

Rev.0 (19/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it