

SUPER

Compact universal head immobilizer



Super head restraints are accessories for spine boards or pick-up stretchers, to be used to increase the degree of patient immobilisation.

Specific features

- The padding of the base increase patient comfort
- Wide strap surfaces make the headblock positioning easier
- Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
- Pillows made of non-absorbent material for an easy sanitation
- Large holes for ear inspection

Technical data

Base dimensions ^(buttonhole excluded)	365 x 260 x 20 mm
Pillow dimensions ^(1 piece)	255 x 85 x 160 mm
Hole for ear inspection	Ø 70 mm
Base weight	360 ± 10 g
Pillow weight	500 ± 10 g
Overall weight	1,46 ± 0,07 kg
Materials	PE, PVC, Nylon

Class I MD compliant with UE Reg. 2017/745

SUPER BLUE YELLOW - COMPACT
UNIV. HEAD IMMOB.
SH00111C



CND Classification	Y060312
Registration number	161554

Rev.0 (19/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it