

755

Universal disposable head immobilizer



Disposable head immobilizer characterized by high versatility at a low cost.

**Specific Characteristics**

- Easy fixing thanks to the integrated adhesive labels
- Soft nape support increases comfort
- Padded adjustable strap for the application of the device

**Technical data**

Dimensions	400x130x10 mm
Materials	PE, PU
Weight	55 g

Class I Med.Device according to Dir. 93/42/CEE

755 - UNIVERSAL DISPOSABLE  
HEAD IMMOBILIZER

**SH00130C**

**CND Classification** Y060312

**Registration number** 194516



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)