

## FXA PRO

### Anatomic universal ankle immobilizer



FXA Pro ankle restraints are accessories for spine boards or pick-up stretchers, to be used to increase the degree of patient immobilisation.

#### Specific features

The wrap-around shape of the support area makes it imperceptible to the patient

One-piece adjustable strap band for device fixing and for ankle immobilization

The band can be removed allowing easy cleaning procedures

Extremely robust and lightweight

#### Technical data

Dimensions	290 x 160 x 80 ± 5 mm
Materials	PE, Nylon
Weight	350 ± 10 g

Class I MD compliant with UE Reg. 2017/745

FXA PRO - ANATOMIC UNIVERSAL  
ANKLE IMMOBILIZER

SH00151A

CND Classification V0899

Registration number 195166



Rev.0 (19/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)