

751

Inflatable disposable head immobilizer



Disposable head immobilizer characterized by high versatility at a low cost.

Specific Characteristics

The transparency limits patient's fear and sense of claustrophobia

Can be fixed with surgical tape

The anti-deflating system integrated in the valve allow an easy adjustment of the inflating pressure making the device extremely adaptable to the patient's conformation
Arrows for an immediate identification of the direction of application

Extremely low storage space

Technical data

Pillow dimensions ^(one piece)	300x200x100 ± 10 mm
Material	PVC
Overall weight	150 ± 10g

Class I Med.Device according to Dir. 93/42/CEE

751 - INFLATABLE DISPOSABLE
HEAD IMMOBILIZER

SH00160A

CND Classification Y060312

Registration number 194499



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it