

## SUPER BLUE HP

### Universal head immobilizer



Super Blue HP head restraints are accessories for spine boards or pick-up stretchers, to be used to increase the degree of patient immobilisation.

#### Specific features

- Wide strap surfaces make the headblock positioning easier
- On the lower side of the base, a large strap element covers the entire surface allowing a more stable fixing on the immobilization device
- Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
- Pillows made of non-absorbent material for an easy sanitation
- Large holes for ear inspection

#### Technical data

Base dimensions	445 x 280 x 15 mm
Pillow dimensions <sup>(1 piece)</sup>	255 x 80 x 155 mm
Hole for ear inspection	∅ 70 mm
Base weight	430 ± 10 g
One pillow weight	500 ± 30 g
Overall weight	1,56 ± 0,3 kg
Materials	PE, PVC, Nylon

Class I MD compliant with UE Reg. 2017/745

SUPER BLUE HP UNIVERSAL HEAD  
IMMOBILIZER  
SH00162C



CND Classification	Y060312
Registration number	477934

Rev.0 (19/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)