

SUPER YELLOW HP

universal head immobilizer



Device for head immobilization intended to be used with Spencer spinal boards or pickup stretchers. It can be used along with other immobilization devices like SED and all Spencer cervical collars.

Specific Characteristics

- Wide strap surfaces make the headblock positioning easier
- On the lower side of the base, a large strap element covers the entire surface allowing a more stable fixing on the immobilization device
- Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
- Pillows made of non-absorbent material for an easy sanitation
- Large holes for ear inspection

Technical data

Base dimensions	445x280x15 mm
Pillow dimensions ^(1 piece)	255x80x155 mm
Hole for ear inspection	ø70 mm
Base weight	430 ± 10g
One pillow weight	500 ± 30g
Overall weight	1,56 ± 0,3 kg
Materials	PE, PVC, Nylon

Class I Med.Device according to Dir. 93/42/CEE

SUPER YELLOW HP FERMACAPO
UNIVERSALE
SH00163C



CND Classification	Y060312
Registration number	801081

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it