

SUPER BLUE HP

Universal head immobilizer



Super Blue HP head restraints are accessories for spine boards or pick-up stretchers, to be used to increase the degree of patient immobilisation.

Specific features

- Wide strap surfaces make the headblock positioning easier
- On the lower side of the base, a large strap element covers the entire surface allowing a more stable fixing on the immobilization device
- Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
- Pillows made of non-absorbent material for an easy sanitation
- Large holes for ear inspection

Technical data

Base dimensions	445 x 280 x 15 mm
Pillow dimensions ^(1 piece)	255 x 80 x 155 mm
Hole for ear inspection	∅ 70 mm
Base weight	700 ± 10 g
One pillow weight	500 ± 30 g
Overall weight	1,83 ± 0,3 kg
Materials	PE, PVC, Nylon

Class I MD compliant with UE Reg. 2017/745

MOD.F011 SUPER BLUE HP
UNIVERSAL HEAD IMMOBILIZER
SH00166C



CND Classification

Y060312

Rev.0 (19/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it