

SPENCER CONTOUR

Anatomic universal head immobilizer



Spencer Contour head restraints are accessories for spine boards or pick-up stretchers, to be used to increase the degree of patient immobilisation.

Specific features

- Openings for ear inspection
- Soft pillows to increase comfort
- Wide strap surfaces make the headblock positioning easier
- Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
- Smooth surfaces for easy sanitation
- Central removable padded area for adult or pediatric use

Technical data

Pillow dimensions ¹	235 x 130 x 160 mm
Base dimensions ²	390 x 260 x 15 mm
Head immobilizer weight	1,5 ± 0,1 kg
Materials	PU, Nylon, PVC

¹One pillow dimensions

² Buttonhole excluded

Class I MD compliant with UE Reg. 2017/745

SPENCER CONTOUR UNIVERSAL
HEAD IMMOBILIZER
YELLOW/BLACK

SH00201A

CND Classification Y060312

Registration number 102237



Rev.0 (19/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it