

## CHIN AND HEAD BELTS



Couple of chin and front spare straps to be used with head immobilizers Spencer non- HP.

### Specific features

Belts with adjustable locking strap for maximum versatility and adaptability

Contact zones realized with soft and high adhesion material to optimize immobilization

### Technical data

Extended belt length <sup>1</sup>	110 cm
Belt width	2 cm
Contact area dimensions	170x35 mm
Weight	90 ± 5g

<sup>1</sup> The useful length varies depending on the head immobilizer in use

CHIN/HEAD BELTS KIT FOR HEAD IMMOBILIZERS NO HP  
SH00219A

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)