

TANGO FIX

Integrated adult/paediatric head immobilizer



Tango Fix head restraints are accessories for spine boards or pick-up stretchers, to be used to increase the degree of patient immobilisation.

Specific features

- The base padding offers high patient comfort
- The central part of the base is removable to allow a better adaptation when the device is used in paediatric configuration
- Base equipped with PVC cover for an effective sanitation
- Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
- The paediatric head immobilizer is obtained from the adult one easily, ensuring rapidity in rescue operation and minimizing storage space inside the ambulance
- The rigid pillows with padded internal side, allow high degree of immobilization without compromising patient comfort
- Adult and paediatric head immobilizer, have large side openings to allow ear inspection
- The pillows are fixed on the base by means of wide strap surfaces

Spare parts

SH00215A	Base for Spencer Contour/Tango fix
SH00219A	Chin/head belts kit for head immobilizers (not HP)
SH00221A	Chin/head belts kit for paediatric head immobilizer

Technical data

Overall dimensions	380 x 265 x 180 ± 10 mm
Paediatric pillow dimensions	190 x 70 x 120 ± 5 mm
Materials	PE, Nylon, TPU
Base weight	450 ± 10 g
Paediatric pillows weight	560 ± 10 g
Adult pillows weight	1060 ± 20 g
Overall weight	2080 ± 30 g

Class I MD compliant with UE Reg. 2017/745

TANGO FIX - ADULT/PAEDIATRIC
HEADIMMOBILIZER

SH00240A

CND Classification



Y060312

Registration number

456414

Rev.0 (19/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it