

TANGO FIX

Integrated adult/pediatric head immobilizer



Device for head immobilization intended to be used with all Spencer spinal boards except pediatric models.

Specific Characteristics

- The base padding offers high patient comfort
- The central part of the base is removable to allow a better adaptation when the device is used in pediatric configuration
- Base equipped with PVC cover for an effective sanitation
- Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
- The pediatric head immobilizer is obtained from the adult one easily, ensuring rapidity in rescue operation and minimizing storage space inside the ambulance
- The rigid pillows with padded internal side, allow high degree of immobilization without compromising patient comfort
- Adult and pediatric head immobilizer, have large side openings to allow ear inspection
- The pillows are fixed on the base by means of wide strap surfaces

Technical data

| | |
|-----------------------------|---------------------|
| Overall dimensions | 380x265x180 ± 10 mm |
| Pediatric pillow dimensions | 190x70x120 ± 5mm |
| Materials | PE, Nylon, TPU |
| Base weight | 450 ± 10 |
| Pediatric pillows weight | 560 ± 10g |
| Adult pillows weight | 1060 ± 20 g |
| Overall weight | 2080 ± 30 g |

Class I Med.Device according to Dir. 93/42/CEE

TANGO FIX - ADULT/PAEDIATRIC HEAD
IMMOBILIZER

SH00240A

CND Classification Y060312

Registration number 456414



Spare parts

| | |
|----------|---|
| SH00215A | Base for spencer contour/tango fix |
| SH00219A | Chin/head belts kit for head immobilizers No Hp |
| SH00221A | Chin/head belts kit pediatric head immobilizer |

Rev.3 (22/10/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it