

## PEDI FIX

### Paediatric head immobilizer



Pedi Fix head restraints are accessories for spine boards or pick-up stretchers, to be used to increase the degree of patient immobilisation.

#### Specific features

- Wide straps allow to place the pillows easier
- On the lower side of the base, a wide strap area covers the entire surface, allowing a more stable fixation of the device
- Quick and easy to apply thanks to the straps on the base
- Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
- The rigid pillows with padded internal side, allow high degree of immobilization without compromising patient comfort
- Large side openings to allow ear inspection

#### Spare parts

- |          |   |
|----------|---|
| SH00221A | Chin/head belts kit paediatric head immobilizer |
|----------|---|

#### Technical data

Pillow dimensions <sup>(one piece)</sup>	190 x 70 x 120 ± 5 mm
Base dimensions	445 x 280 x 15 ± 5 mm
Materials	PE, Nylon, TPU
Pillow wight <sup>(one piece)</sup>	560 ± 10 g
Base weight	430 ± 10 g

Class I MD compliant with UE Reg. 2017/745

PEDI FIX PAEDIATRIC  
HEADIMMOBILISER

SH00260A

CND Classification Y060312

Registration number 1069547



Rev.0 (19/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)