

PEDI GO

Paediatric head immobilizer



Pedi Go head restraints are accessories for spine boards or pick-up stretchers, to be used to increase the degree of patient immobilisation.

Specific features

Quick and Easy to apply thanks to the strap on the bottom
 Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
 The rigid pillows with padded internal side, allow high degree of immobilization without compromising patient comfort
 Large side openings to allow ear inspection

Spare parts

SH00219A	Chin/head belts kit for head immobilizers (not HP)
SH00221A	Chin/head belts kit paediatric head immobilizer

Technical data

Pillow dimensions ^(one piece)	190 x 70 x 120 ± 5 mm
Materials	PE, Nylon, TPU
Paediatric pillows weight	560 ± 10 g

Class I MD compliant with UE Reg. 2017/745

PEDI GO PAED. HEAD IMM. CUSH. FOR
 BABY GO/PEDI LOC

SH00262A

CND Classification



Y060312

Registration number

425600

Rev.0 (19/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it