

## SUPER SX

### Head immobilizer for pick-up stretcher



Super SX head restraints are accessories for pick-up stretchers, to be used to increase the degree of patient immobilisation.

#### Specific features

- The padding of the base increase patient comfort
- Wide strap surfaces make the headblock positioning easier
- Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
- Pillows made of non-absorbent material for an easy sanitation
- Large holes for ear inspection
- The base split in two pieces, allows to use the head immobilizer without the need to remove it during the application of the pick-up stretcher

#### Technical data

Base dimensions <sup>(1 half)</sup>	200 x 190 x 20 mm
Pillow dimensions <sup>(1 piece)</sup>	255 x 85 x 160 mm
Hole for ear inspection	∅ 70 mm
Base weight	360 ± 10 g
Pillow weight	500 ± 10 g
Overall weight	1,46 ± 0,07 kg
Materials	PE, PVC, Nylon

Class I MD compliant with UE Reg. 2017/745

SPENCER SX-HEAD IMMOBILIZER  
FOR PICK UP STRETCHER

SH00500

CND Classification



Y060312

Rev.0 (19/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)