

## Pro Skid-E AIR

### Evacuation/transport chair



Pro Skid – E AIR is an evacuation chair ideal for transfers on stairs thanks to the front and rear handles which allow lifting, and slides that allow easy movement on stairs. 4 wheels allow movement on flat surfaces. This particular version is characterized by reduced width, making it suitable for use on aircrafts

#### Specific features

- PVC Sheet easy to clean
- Slide allow easy movement and control going downstairs
- Telescopic front handles and rear handles make going upstairs easier
- Footrest
- Front wheels pivoting with brake
- Large diameter rear wheels increase sliding action
- Headrest sheet is movable making the device suitable for patient of different heights

#### Reference standard

N.A.

#### Standard equipment

Couple of belts

#### Optionals

- |          |                               |
|----------|-------------------------------|
| SK21000E | Skid carry bag                |
| SK14000B | WALL BRACKET FOR SKID/PROSKID |

#### Technical data

Height opened	1580 mm
Height folded	1060 mm
Length with handles closed	1110 mm
Length with handles opened	1340 mm
Width	410 mm
Thickness folded	355 mm
Weight	14 kg
Loading capacity	150 kg
Materials	Al, PVC, PE,PU
Front wheels	ø100 mm
Rear wheels	ø200 mm

Dimensions and weight are subjected to tolerances of ±10mm and ± 1kg

Class I Med.Device according to Dir. 93/42/CEE

PROSKID-E AIR YELLOW  
FRAME, BLACK COVER



SK10002E

CND Classification

V0899

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)