

Pro Skid-E MAX

Evacuation/transport chair



Pro Skid – E MAX is an evacuation chair ideal for transfers on stairs thanks to the front and rear handles which allow lifting, and slides that allow easy movement on stairs. 4 wheels allow movement on flat surfaces. This particular version is characterized by increased loading capacity.

Specific Characteristics

- PVC Sheet easy to clean
- Slide allow easy movement and control going down and upstairs
- Telescopic front handles make handling easier
- Footrest
- Front wheels pivoting with brake
- Large diameter rear wheels increase sliding action
- Headrest sheet is movable making the device suitable for patient of different heights

Standard equipment

ST00427A BLACK PAIR OF BELTS

Accessories

SK21000E Skid carry bag
SK14000B WALL BRACKET FOR SKID/PROSKID

Technical data

Height opened	1580 mm
Height folded	1060 mm
Length with handles closed	1110 mm
Length with handles opened	1340 mm
Width	550 mm
Thickness folded	355 mm
Weight	15 kg
Loading capacity	250 kg
Materials	Al, PVC, PE, PU
Front wheels	Ø100 mm
Rear wheels	Ø 200mm

Dimensions and weight are subjected to tolerances of ±10mm and ± 1kg

Class I Med.Device according to Regulation 2017/745/UE

PROSKID-E MAX YELLOW
FRAME, BLACK SHEET
SK10008E



CND Classification

V0899

Rev.1 (22/09/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it