

## Pro Skid-E

### Evacuation/transport chair



Pro Skid – E is an evacuation chair ideal for transfers on stairs thanks to the front and rear handles which allow lifting, and slides that allow easy movement on stairs. 4 wheels allow movement on flat surfaces.

#### Specific Characteristics

PVC Sheet easy to clean
Slide allow easy movement and control going down and upstairs
Telescopic front handles make handling easier
Footrest
Front wheels pivoting with brake
Large diameter rear wheels increase sliding action
Headrest sheet is movable making the device suitable for patient of different heights

#### Standard equipment

ST00427A	BLACK PAIR OF BELTS
----------	---------------------

#### Accessories

SK21000E	Skid carry bag
SK14000B	WALL BRACKET FOR SKID/PROSKID

#### Technical data

Height opened	1580 mm
Height folded	1060 mm
Length with handles closed	1110 mm
Length with handles opened	1340 mm
Width	550 mm
Thickness folded	355 mm
Weight	15 kg
Materials	Al, PVC, PE,PU
Front wheels	ø100 mm
Rear wheels	ø 200mm

Dimensions and weight are subjected to tolerances of ±10mm and ± 1kg

#### Class I Med.Device according to Dir. 93/42/CEE

PROSKID-E SILVER FRAME,  
BLACK SHEET  
SK10201E



CND Classification	V0899
Registration number	129682

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)