

## SKID WALL BRACKET



Wall bracket to fix on the wall an evacuation chair of SKID series, compatible with model SKID-E / SKID-E PRO / SKID-E READY , not compatible with SKID-OK.

### Specific Characteristics

Simple to fix the bracket on the wall (dowels included 8x50mm)

Simple to fix the chair on the bracket

### Technical data

Dimension	365x185x170mm ±5mm
Material	Steel
Weight	0,4kg±0,1kg

WALL BRACKET FOR SKID/PROSKID  
**SK14000B**

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)