

## SKID OK WALL BRACKET



Wall fastener for SKID OK and SKID OK B evacuation chairs.

### Specific features

Easy to install (dowels 8x50mm included)

The placing of the chair on the fastener is very simple

### Technical data

Dimensions	272x166x164±5mm
Material	Steel
Weight	0,4 ± 0,1kg

WALL FASTENER FOR SKID OK

**SK14001B**

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)