

Skid-E

Evacuation chair



Skid – E is an evacuation chair ideal for transfers on stairs thanks to the integrated slides that allow easy movement on stairs. 4 wheels allow movement on flat surfaces.

Specific Features

PVC Sheet easy to clean
Slide allow easy movement and control going down and upstairs
Footrest
Front wheels pivoting with brake
Large diameter rear wheels increase sliding action
Headrest sheet is movable making the device suitable for patient of different heights

Standard equipment

ST00427A	BLACK PAIR OF BELTS
----------	---------------------

Accessories

SK21000E	Skid carry bag
SK14000B	WALL BRACKET FOR SKID/PROSKID

Technical data

Height opened	1580 mm
Height folded	1060 mm
Width	530 mm
Length	900 mm
Thickness folded	330 mm
Weight	12,7 kg
Loading capacity	150 kg
Materials	Al, PVC, PE,PU
Front wheels	ø100 mm
Rear wheels	ø 200 mm

Dimensions and weight are subjected to tolerances of ±10mm and ± 1kg

Class I Med.Device according to Dir. 93/42/CEE

SKID-E SILVER FRAME BLACK

SEAT

SK20001E



CND Classification V0899

Registration number 129686

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it